



Introduction

For most therapeutic glycoproteins the glycosylation patterns greatly influence clinical performance of the drug product, particularly its *in vivo* safety and efficacy profile.¹ In biological tissues glycosylation patterns can also correlate with the state of health or disease of the individual.² Given this, there is increasing interest in accurately characterizing glycosylation, for example monitoring glycosylation patterns of biopharmaceutical therapeutics throughout the product lifecycle as well as in glycan biomarker discovery for medical diagnostics.

Robust analytical strategies are required to meet the challenge of accurately and reliably characterizing glycosylation. There has been significant progress made in glycan analysis and the availability of commercial kits which contain the necessary reagents for release and labelling of monosaccharides, sialic acids, and *N*- and *O*-glycans have made it easier for laboratories to adopt technologies for glycan analysis. However, even with the advancement in glycan characterisation tools, multiple inter-laboratory studies have shown that there is still a lack in consistency of the data produced during glycan analysis.³ These problems highlight the existing need for well-characterised glycan reference standards.

Glycan Standards

A key component in a well-designed analytical strategy is the inclusion of standards. Table 1 shows which standards can be used for best practice during the analysis of sialic acids, monosaccharides, *N*-glycans and/or *O*-glycans. These fall into the following categories;

1. System suitability standards enable an analyst to test the holistic functionality of an analytical system (e.g. chromatographic, mass spectroscopic and/or CE) and evaluate whether it is adequate for its intended use.

2. Process standards or process controls are used to verify that part of or an entire process has worked correctly. There are four main categories for processes standards in glycoanalysis: release, labelling, release followed by labelling and exoglycosidase sequencing

3. Reference standards allow for characterisation by comparison. This can be accomplished by the direct comparison of the chromatographic or electrophoretic retention time of an unknown to that of a standard whose structure have been fully characterised. Additionally, primary assignment of unknown structures can be accomplished by comparison of their GU (Glucose unit) values (obtained using glucose homopolymer (GHP) standard) with glycans whose GU values are in databases or in the literature.

4. Quantitative glycan standards are used to determine the absolute amount of an analyte in a sample. They can also be used to quantify the efficiency of a process.

System Suitability Standards and Controls		Process positive control for release/labelling analysis				Release Process positive control				Labelling Process positive control				MS System suitability				LC/HPLC System suitability				GC Calibration				WAX System Suitability				CE System Suitability				Structure Identification by retention time or m/z matching				Quantification				Exoglycosidase control			
		SA	MONO	N	O	SA	MONO	N	O	SA	MONO	N	O	N	O	SA	MONO	N	O	N	O	N	O	N	O	N	O	SA	MONO	N	O	SA	MONO	N	O	SA	MONO	N	O						
Glycoprotein/Glycopeptides	Ludger Code																																												
	GCP-IGG	IgG Glycoprotein																																											
	GCP-FET	Fetuin Glycoprotein																																											
Monosaccharide and Sialic Acid standards	BQ-GPEP-A2G2S2	GPEP-A2G2S2																																											
	CM-SRP	SRP Sialic Acid Reference Panel																																											
	CM-NEU-AC	Neu5Ac																																											
	CM-NEU-GC	Neu5Gc																																											
	CM-NEU,9AC2	Neu5,9Ac2																																											
	CM-MONOMIX	MonoMix																																											
Unlabelled N-Glycans	CM-XYL	Xylose																																											
	CN-x	Bi, Tri and Tetra-antennary N-glycans																																											
	CN-Man-x	High Mannose N-glycans																																											
	BQ-CNTOTRIOSE	Chitotriose																																											
	BQ-CN-MAN8	Man8																																											
	CLBN-IGG	IgG N-glycan library																																											
	SA-MAB4	Mab 4 glycan ref panel																																											
2-AB labelled glycans	CLBN-FETUIN	Fetuin N-glycan library																																											
	CLBO-FETUIN	Fetuin O-glycan library																																											
	CAB-GHP	2-AB labelled GHP																																											
2-AA labelled glycans	BQ-CAB-CHI	2-AB labelled Chitotriose																																											
	CAB-IGG	2-AB labelled IgG N-glycan library																																											
	CAB-x	Bi, Tri and Tetra-antennary N-glycans																																											
	CAB-Man-x	High Mannose N-glycans																																											
	CAB-C-x	O-glycans																																											
	CAB-AlphaGal	Alpha-Gal standard																																											
APTS labelled glycans	CAA-GHP	2-AA labelled GHP																																											
	BQ-CAA-CHI	2-AA labelled Chitotriose																																											
	CAA-x	Bi, Tri and Tetra-antennary N-glycans																																											
	CAA-Man-x	High Mannose N-glycans																																											
PROC labelled glycans	CAA-AlphaGal	Alpha-Gal standard																																											
	CAPT5-IGG	APTS labelled IgG N-glycan library																																											
Permethylyated Glycan Standards	CAPT5-x	N-glycans																																											
	CPRC-GHP	PROC labelled GHP																																											
	CPRC-IGG	Proc labelled IgG N-glycan library																																											
	CPRC-x	Bi and Tri-antennary N-glycans																																											
Permethylyated Glycan Standards	CPRC-Man-x	High Mannose N-glycans																																											
	CFM-IGG	Permethylyated IgG N-glycans																																											
	CFM-C13-IGG	Permethylyated (13C) IgG N-glycans																																											

Table 1: Standards Used for Glycan Analysis






Quantitative Analysis: Case Studies for Glycan Standards

Here we show how various glycan standards work in concert to provide confidence in results. As case studies, we will illustrate how we use each type of standard to support reliable and consistent sialic acid analysis and monosaccharide analysis. These methods are used to quantify the total mass of specific monosaccharides within the glycan pool. Both methods are required to satisfy regulatory requirements for biopharmaceutical drug characterisation (e.g. ICH guideline Q6B).

1) Sialic Acid Analysis

Sialylation is integral to the structure and function of many glycoproteins and is a glycosylation critical quality attribute (GCQA) for biopharmaceuticals. Sialic acid analysis provides data for both the abundance and the type of sialylation (including O-acetylation).

Workflow:



Sialic Acid	Amount of Sialic Acids on Protein (nmoles/mg protein)	Amount of Neu5Gc and Neu5Ac (nmoles)	Average Relative Percent (%) of Sialic Acids from Peak Areas
Neu5Gc	4.97	0.25	0.95
Neu5Ac	307.44	15.37	94.17
Neu5,9Ac2	-	-	4.88

✓ **System suitability Standard and Reference Standard**
Sialic acid reference panel (SRP) - contains a mixture of sialic acids found in humans and animals.
Acceptance criteria: the HPLC profiles at the start and end of the sample set should overlap with minimal drift in retention time (e.g. ± 0.1 min.).

✓ **Quantitative Standards**
Neu5Gc and Neu5Ac - dispensed to approximately 1 nmole. Preparation of serial dilutions provide calibration curves.
Acceptance criteria: the calibration curves should give R^2 values of >0.99



✓ **Quantitative Process Standard**
GPEP-A2G2S2 glycopeptide - contains a biantennary *N*-linked glycan terminating in two sialic acids and has been quantified using qNMR.
Acceptance criteria: Neu5Ac in the range of 5.6 to 8.4 nmol (which is the amount determined by quantitative NMR $\pm 20\%$)

✓ **Reliable Data for Erythropoietin (EPO):** EPO is a highly glycosylated protein and its high level of sialylation and accompanying acetylation has a significant effect on its therapeutic properties.⁴ The DMB labeled sialic acids from an EPO glycoprotein were analysed and the relative levels of the *N*-acetyl, *N*-glycolyl and *O*-acetyl sialic acids were calculated as a ratio of 94:1:5 respectively.



2) Monosaccharide Analysis

Monosaccharide analysis provides absolute or relative quantitation of the neutral (i.e. non-anionic) monosaccharides and information relating to the types of *N*- and/or *O*-glycans present on a glycoprotein.

Workflow: Monosaccharides are released from the glycoprotein by acid hydrolysis using 2M trifluoroacetic acid (TFA) for quantitation of mannose (Man), galactose (Gal), glucose (Glc), fucose (Fuc) or 6M hydrochloric acid (HCl) for quantitation of glucosamine (GlcN) and galactosamine (GalN)



Monosaccharide (HCL & TFA Data)	Average (nmol)	%CV
GlcN	15.46	1.0
GlcN	9.80	3.4



HCL & TFA Release Data	nmol/mg protein	%CV	# of monosaccharides per N-glycan site
GlcN	58	2.1	4.4
GalN	3	2.0	0.2
Gal	24	2.4	1.8
Man	41	1.6	3.1
Glc	5	9.2	0.4
Fuc	6	1.3	1.2

✓ **Reference Standard, Quantitative Standard, System Suitability Standard**
Monosaccharide mix (MonoMix) standard - contains GlcN, GalN, Gal, Man, Glc and Fuc dispensed to 10 nmole each. Preparation of serial dilutions provide calibration curves.
Acceptance criteria for quantitation: calibration curves should give R^2 values of >0.9
System suitability acceptance criteria: HPLC retention times of the monosaccharides should have less than 0.1 min difference when run at the start, middle and end of the analysis.

✓ **Quantitative Process Standards**
GPEP-A2G2S2 glycopeptide or **Man-8** quantified glycan standards
Acceptance criteria for GPEP-A2G2S2: GlcN is in the range of 8.38 to 13.96 nmol (TFA) and 11.17 to 16.75 nmol (HCl).

✓ **Reliable Data for Monoclonal Antibody (mAb):** Glycosylation is a common post-translational modification in mAbs and has a critical role in antibody effector function.⁵ The monosaccharide analysis of a mAb gave the following;

- From the data for GlcN and GalN (from the HCl release):
- majority of the *N*-glycans are biantennary (with 4 GlcN)
- low percentage of *O*-glycans (low, but real GalN value)

- From the data for Gal, Man, Glc & Fuc (from TFA release)
- Not all *N*-glycans have 2 Gal (1.8 per site).
- Complex *N*-glycans have 3 Man; indicates low % of high mannose structures present.
- The majority of the *N*-glycans are core fucosylated (1.2 per site)
- Glc is present as a background contaminant (detected in the negative controls)

Acknowledgements and References

The human EPO was a kind gift from Raquel Montesino and Antonio Vallin from the Center for Genetic Engineering and Biotechnology, Cuba.

- P. Zhang, S. Woen, T. Wang, B. Liaw, S. Zhao, C. Chen, Y. Yang, Z. Song, M.R. Wormald, C. Yu, P.M. Rudd. *Drug Discov Today*, 2016, **21**,740-765.
- G. Lauc, M. Pezer, I. Rudan, H. Campbell. *Biochim Biophys Acta - Gen Subj*, 2016, **1860**, 1574-1582; B. Adamczyk, T. Tharmalingam, P.M. Rudd. *Biochim Biophys Acta - Gen Subj*, 2012, **1820**,1347-1353; A. Varki. *Glycobiology*, 2017, **27**,3-49.
- S. Thobhani, C.T. Yuen, M.J.A. Bailey, C. Jones. *Glycobiology*, 2009, **19**,201; M. L. A. D. Leoz, D. L. Dweyer, S. E. Stein, NISTIR 8186 NIST Interlaboratory Study on the Glycosylation of NISTmAb, a Monoclonal Antibody Reference Material June 2015 to February 2016, 2016.
- M.N. Fukuda, H. Sasaki, L. Lopez, M. Fukuda, M. Blood, 1989, **73**, 84-89.
- K. Zheng, C. Bantog, R. Bayer, *mAbs*, 2011, **3**, 568-576

Contact for more information: If you'd like a copy or want to know more please email me (jenifer.hendel@ludger.com) or connect with me on LinkedIn.

