

# LudgerTag V-Tag for N-glycan profiling and identification



Ludger





# Highlights of the V-Tag System

## V-Tag labeling

*Allows glycan identification and quantitation using (U)HPLC.  
Provides data comparable to gold-standard glycoprofiling method  
based on 2-AB. V-Tag dye is 10 times more sensitive than 2-AB.*

### Sample Needed

*V-Tag labeling is suitable for  
10 $\mu$ g to 200 $\mu$ g of sample in a  
IgG type sample*

### Validated for GMP Labs

*Validated to ICH Q2(R1) standards  
using IgG and Fetuin (Highly  
sialylated type sample)*



### Quick and Easy

*PNGase F, labeling and clean up  
is completed within 2 hours for  
30 IgG type samples.*

### Applications

*This technology can be used for  
biopharma studies (e.g. IgG) and  
different biological samples  
(e.g. plasma)*


### Exoglycosidase digestion

*V-Tag labeled glycans are compatible with exoglycosidase  
digestion, including the removal of core alpha (1-6) fucose by  
Fucosidase from Bovine kidney (Fuc-BK)*

# Anatomy of V-Tag

The molecular aspects that make V-Tag work

**Fluorescent** group for  
detection in UHPLC  
 $\lambda_{\text{ex}} = 250\text{nm}$ ,  $\lambda_{\text{em}} = 360\text{nm}$



**Amine reactive succinimidyl ester**  
to react with the amino glycans

**Simple, non-reactive alkyl chain** to  
link the fluorescent moiety to the reactive  
succinimidyl ester

# Components of the LT-VTAG-C30 Kit


## 1. Glycan release



*Glycan release kit*

**LudgerZyme™ rPNGase F release kit**  
LZ-rPNGaseF-30

## 2. Labeling




*Reaction solvent*

**DMSO**  
LT-DMSO-02

*Amine reactive fluorescent labeling reagent*

**V-Tag Labeling Dye**  
LT-VTAG-02

## 3. Clean up




*HILIC resin cartridge*

**LudgerClean™ VSPE Cartridges**  
LC-VSPE-30

*Solvents involved in HILIC clean-up*

**ammonium formate Buffer, pH 4.4**  
LC-N-BUFFX40-30  
**ammonium acetate Buffer**  
LC-AA-BUFF-30

# Workflow - LudgerTag V-Tag for N-glycan profiling and identification



Glycoprotein

## Glycan release

1. Denaturation step (10min, 100°C)
  2. PNGase F (10min, 37°C)
- (Up to 3 hours for more complex samples)

N-glycans

## Labeling

V-Tag Dye  
(15min, 37°C)

## Clean up

SPE clean up (1h)

V-Tag labeled glycans

*Labels the amino group of the released glycans*


Symbol for sugar

- Glc
- GlcNAc
- ★ NeuNAc
- ◇ Gal
- ◆ GalNAc
- ◇ Fuc (deoxy galactose)
- Man

Fluorescent label

—PROC Procainamide

Linkage position




Linkage type

- β-linkage
- ..... α-linkage
- ~ unknown β-linkage
- ~ unknown α-linkage

# Sample preparation within a day-N-glycans

## Timing for processing of 30 samples





# Typical data



(U)HPLC profile


ESI-MS profile



# Repeatability

V-Tag System Has Been Validated to ICH Q2(R1) Level

## IgG type samples




Overlaid IgG V-Tag labeled glycan profiles of nine UPLC analyses at 3 different concentrations (25µg, 50µg and 200µg)

Peak Number	Relative Area %		
	Av.	Std. Dev	CV
1	19.07	0.58	3.05
2	5.45	0.14	2.58
3	19.52	0.33	1.67
4	8.51	0.23	2.76
5	4.99	0.05	0.97
6	1.06	0.09	8.04
7	14.37	0.17	1.20
8	1.68	0.11	6.44
9	2.80	0.11	3.82
10	2.29	0.18	7.77
11	10.72	0.23	2.16
12	2.39	0.14	5.66
13	1.60	0.12	7.35
14	2.72	0.14	5.28
15	2.83	0.12	4.41

Average relative % area, SD and CVs for V-Tag labelled IgG glycans

Validation study shows repeatability for relative % areas

## Highly sialylated type sample



Overlaid Fetuin V-Tag labeled glycan profiles of nine UPLC analyses at 3 different concentrations (25µg, 50µg and 200µg)



Peak Number	Relative Area %		
	Av.	Std. Dev	CV
1	5.82	0.12	2.03
2	3.28	0.09	2.73
3	2.95	0.15	5.08
4	33.48	0.39	1.16
5	23.50	0.27	1.13
6	20.14	0.52	2.56
7	6.68	0.11	1.63
8	1.55	0.05	3.51
9	1.77	0.03	1.63
10	0.83	0.04	4.20

Average relative % area, SD and CVs for V-Tag labeled Fetuin glycans



# Intermediate Precision

## IgG type sample




Overlaid IgG V-Tag labeled glycans profiles from day 1 (A) and day 2 (B) both run in triplicate (90µg of glycoprotein used)

The variation between samples when processed on two separated days

Peak number	Relative Area %					
	Day 1 average	Day 2 average	Day 1 log Av	Day 2 log Av	Ratio Av	% Variation
1	18.38	18.03	1.26	1.26	1.02	1.94
2	5.46	5.84	0.74	0.77	0.93	6.58
3	18.79	18.71	1.27	1.27	1.00	0.44
4	8.31	8.24	0.92	0.92	1.01	0.85
5	4.84	4.89	0.68	0.69	0.99	1.02
6	1.06	1.22	0.02	0.09	0.86	13.58
7	14.07	13.58	1.15	1.13	1.04	3.60
8	1.56	1.77	0.19	0.25	0.88	11.60
9	2.93	3.07	0.47	0.49	0.95	4.79
10	2.97	3.34	0.47	0.52	0.89	11.20
11	10.96	10.75	1.04	1.03	1.02	1.89
12	2.47	2.44	0.39	0.39	1.01	1.16
13	1.99	2.14	0.30	0.33	0.93	7.01
14	3.06	2.90	0.49	0.46	1.06	5.55
15	3.17	3.07	0.50	0.49	1.03	3.08



Average relative % area, SD and CVs for V-Tag labeled IgG glycans from day 1 and 2



IgG V-Tag average relative area % from day 1 compared to day 2. LT-VTAG-C30 9

# Intermediate Precision

## Highly sialylated type sample




Overlaid Fetuin V-Tag labeled glycans profiles from day 1 (A) and day 2 (B) both run in triplicate (90µg of glycoprotein used)

The variation between samples when processed on two separated days

Peak number	Relative Area %					
	Day 1 average	Day 2 average	Day 1 log Av	Day 2 log Av	Ratio Av	% Variation
1	5.72	5.87	0.76	0.77	0.97	2.58
2	3.31	3.31	0.52	0.52	1.00	0.01
3	2.89	2.86	0.46	0.46	1.01	1.17
4	32.98	33.43	1.52	1.52	0.99	1.34
5	23.34	23.56	1.37	1.37	0.99	0.91
6	20.67	20.28	1.32	1.31	1.02	1.91
7	6.75	6.54	0.83	0.82	1.03	3.20
8	1.63	1.57	0.21	0.20	1.04	3.90
9	1.85	1.73	0.27	0.24	1.07	7.22
10	0.85	0.85	-0.07	-0.07	1.00	0.23

Average relative % area, SD and CVs for V-Tag labeled Fetuin glycans from day 1 and 2




Fetuin V-Tag average relative area % from day 1 compared to day 2.

# Linearity and Working range


## V-Tag System Has Been Validated to ICH Q2(R1) Level

### IgG type samples




Overlay plot of one replicate for each amount of IgG taken through labeling (10µg, 25µg, 50µg and 200µg)

Stacked plot of one replicate for each amount of IgG taken through labeling.




IgG area counts for peaks 1, 7, 11 and 15 versus the starting amount of IgG for each labeling.

### Highly sialylated type sample



Overlay plot of one replicate for each amount of Fetuin taken through labeling (10µg, 25µg, 50µg and 200µg)

Stacked plot of one replicate for each amount of Fetuin taken through labeling.





Fetuin area counts for peaks 1, 4, 5, 6 and 7 versus the starting amount of Fetuin for each labeling.

# Comparability of V-Tag and 2-AB

IgG type sample

## Relative area comparison



A) UHPLC chromatogram of V-Tag labeled glycans from IgG after using LT-VTAG-30.

B) UHPLC chromatogram of 2-AB labeled glycans from IgG after overnight PNGase F release, 1 hour labeling with 2AB and LudgerClean T1 purification.



**Comparable in quantitation and in reliability, using a HILIC column**



# Comparability of V-Tag and 2-AB

Highly sialylated type sample

## Relative area comparison




A) UHPLC chromatogram of V-Tag labeled glycans from Fetuin after using LT-VTAG-30.

B) UHPLC chromatogram of 2-AB labeled glycans from Fetuin after overnight PNGase F release, 1 hour labeling with 2AB and LudgerClean T1 purification.

**Comparable in quantitation and in reliability, using a HILIC column**

# Exoglycosidase digestion of V-Tag labeled N-glycans



UHPLC chromatogram of V-Tag labeled glycans from IgG after exoglycosidase digestion

**Exoglycosidase digestion of V-Tag labeled N-glycans is possible, including the removal of core alpha (1-6) fucose by Fucosidase from Bovine kidney (Fuc-BK)**

*Und: undigested; Sial: Sialidase from Arthrobacter ureafaciens: specific for  $\alpha$ 2-3, 6, 8, 9 sialic acids [Ludger, E-S001]; Fuc (BK): fucosidase from Bovine kidney, specific for  $\alpha$ 1-6 $\rightarrow$ 2 fucose [Sigma, F-5884]; Fuc (EM): fucosidase from Elizabethkingia miricola: specific for  $\alpha$ 1-6 fucose [E-F006]; Gal:  $\beta$ -galactosidase from Streptococcus pneumoniae: specific for  $\beta$ 1-4 galactose [Ludger, E-BG-07]*

# Next Steps ...

If you have a question



CLICK  
to contact  
Conchi

**Dr Concepcion Badia Tortosa**  
*Senior Scientist*  
[conchi.badia@ludger.com](mailto:conchi.badia@ludger.com)

Request a quotation



CLICK  
to contact  
Sales

**Sales Team**  
[info@ludger.com](mailto:info@ludger.com)